

# Accelerated Stability Assessment Program (ASAP) for the Rapid Determination of Shelf Life for Oligonucleotides

A. K. Waterman<sup>a</sup>, K. Flavie<sup>a</sup>, C. Shende<sup>a</sup>, J. McCabe<sup>a</sup>, C. Reidy<sup>b</sup>, M. Araya<sup>b</sup>, B. Addepli<sup>b</sup>, K. Hall<sup>b</sup>, and K. C. Waterman<sup>a</sup>

<sup>a</sup>FreeThink Technologies, Inc. and <sup>b</sup>Waters Corporation

## Introduction

- Oligonucleotide therapeutics including antisense oligonucleotides (ASOs), small interfering RNAs (siRNAs), and aptamers are generally manufactured as lyophiles or aqueous solutions in vials or pre-filled syringes.
- Predominant degradation pathways include hydrolysis and oxidation.
  - Sensitivity to temperature, humidity, and oxygen has been demonstrated.
- Accelerated Stability Assessment Program (ASAP) case study performed using a lipid-conjugated ASO reference standard provided by Waters Corporation.

5' d 5-Pal-\*<sup>-</sup>-MOE-MeC-\*<sup>-</sup>-MOE-G-\*<sup>-</sup>-MOE-MeC-\*<sup>-</sup>C\*G\*A\*T\*A\*A\*G\*G\*T\*A\*-MOE-MeC-\*<sup>-</sup>-MOE-A-\*<sup>-</sup>-MOE-MeC 3'

\* = phosphorothioate bond (sulfur replacing oxygen in phosphate linkage)


MOE = methoxyethyl 2' modification, MeC = 5-methyl cytidine residue

5' palmitate modification

## Analytical Method

- Loss of active (% purity) monitored using hydrophilic interaction liquid chromatography (HILIC).
  - Hydrophilic (amide) stationary phase packed into MaxPeak™ High Performance Surfaces (HPS): reduces nonspecific binding, eliminates need for column passivation, and increases reproducibility.
  - Oligonucleotides strongly retained at high organic conditions, with retention influenced by strand length, column temperature, and mobile phase ionic strength/aqueous content.


Mobile Phase	A: 90:10 ACN:100 mM NH <sub>4</sub> OAc B: 10:90 ACN:100 mM NH <sub>4</sub> OAc		
Gradient	Time (min)	% A	% B
	0	90	10
	0.67	90	10
	4.67	22	78
	5.00	22	78
	5.33	90	10
13.33	90	10	
Flow Rate	0.3 mL/min		
Column Temp.	75°C		
Injection Vol.	10 µL		



## ASAP Study Design


- Experimental design aims to bracket specification limit (10% loss of active).
- Samples exposed, open, to temperature/relative humidity conditions.

Temp (°C)	Saturated Salt	% RH	Time (d)
Control	NA	NA	0 (n = 3)
25	LiCl	11	12
30	CsF	3	13, 21
35	KF	25	6, 21
40	LiCl	11	4, 9, 21
45	KF	21	2, 5, 16
50	LiCl	11	2, 5, 15
50	NaI	29	1



## Results

- Non-linear behavior observed (two-step process).
- Isoconversion times (time to lose 10% purity) determined using bilinear fitting.



- Isoconversion rates (k) well fit by moisture-modified Arrhenius equation:



$$\ln(k) = \ln A - \frac{E_a}{RT} + B(RH)$$

- In A = intercept, E<sub>a</sub> = activation energy [temp sensitivity], R = gas constant, T = temp in K, B = humidity sensitivity, RH = equilibrium humidity

lnA	E <sub>a</sub> (kcal/mol)	B	R <sup>2</sup>	Q <sup>2</sup>
27.5 ± 3.7	17.2 ± 2.2	0.066 ± 0.011	0.982	0.945

## ASAP Model Predictions

- ASAP shelf life modeling (using ASAPprime®) to 95% probability of passing:
  - Less than six months at 5°C, even very dry (with molecular sieves).
  - More than nine years at -20°C.



## Conclusions

- Three-week ASAP study conducted on lyophilized lipid-conjugated ASO.
  - HILIC method provides effective separation.
- Loss of purity well fit by moisture-modified Arrhenius equation using isoconversion.
  - Lipid-conjugated ASO requires frozen storage (-20°C) for long-term stability.
- Stability of RNA-based molecules is effectively modeled using ASAP.

## Contact

alisa.waterman@freethinktech.com



Scan here to access this poster and more

