

Excipient Compatibility as Predicted by ASAP

A Case Study

Getting To Know Faster

Outline

1. Introduction
2. Overview of Traditional Excipient Compatibility Studies
3. Accelerated Stability Assessment Program (ASAP)
4. Case Study
5. Alternate Analytical Technique
6. Comparison of Analytical Techniques
7. Conclusion
8. Acknowledgement

Introduction

- ▶ Provide a rationale for excipient selection
- ▶ ICH Q8(R2): Pharmaceutical Development
 - ▶ The compatibility of the drug substance with excipients listed should be evaluated
- ▶ Impact on stability of API
 - ▶ Risk management
- ▶ Regulatory Requirement
 - ▶ QOS: What evidence supports compatibility between the excipients and the drug substance?

Traditional Excipient Compatibility Studies

- ▶ Binary mixture of active pharmaceutical ingredient (API) and excipient prepared in a 1:1 ratio
- ▶ Excipients showing incompatibility at the 1:1 ratio tested at a ratio more typical of the final formulation
- ▶ Two vial types (duplicate vials for each) prepared for each drug-excipient blend
 - ▶ Assay vial (pre-weighed for moisture)
 - ▶ Composite vial (for additional tests)

Traditional Excipient Compatibility Study

- ▶ Three controlled conditions
- ▶ 200 vials prepared for 7 excipients studied open dish

One Excipient = 25 Vials	<u>T=0</u>			<u>T=2 weeks</u>			<u>T=4weeks</u>			<u>T=8 weeks</u>		
<u>Test Conditions</u>	<u>Assay</u>	<u>Dupl.</u>	<u>Comp.</u>	<u>Assay</u>	<u>Dupl.</u>	<u>Comp.</u>	<u>Assay</u>	<u>Dupl.</u>	<u>Comp.</u>	<u>Assay</u>	<u>Dupl.</u>	<u>Comp.</u>
Accelerated - 40°C / 75% RH	X	X	X	X	X	X	X	X	X	X	X	X
CRT - 25°C / 60% RH	X	X	X	X	X	X	X	X	X	X	X	X
Light	2 X	X	X	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Total Vial Counts =	4	3	3	2	2	2	2	2	2	2	2	2
Grand Totals =	10			6			6			6		

Traditional Excipient Compatibility Studies

- ▶ The analytical testing includes:
 - ▶ Assay
 - ▶ Impurities
 - ▶ Description
 - ▶ Loss on Drying (LOD)
- ▶ Results compared to a control (API) across all time points and conditions to look for any trending
- ▶ Any unexplained or unexpected loss in potency, or increase in impurities over time mandated further investigations
 - ▶ Excipient ratio screening, at the levels typical of the final formulation.

Traditional Excipient Compatibility Studies

- ▶ Limitations
 - ▶ 1:1 Binary blends not representative of formulation
 - ▶ May not pick up on incompatibilities
 - ▶ Manufacturing
 - ▶ Processing
 - ▶ Tableting
 - ▶ Sample Load
 - ▶ Resources
 - ▶ Increased time-line
 - ▶ Inconclusive results
 - ▶ Not an indicator of formulation stability

Accelerated Stability Assessment Plan: Excipient Compatibility Study

- ▶ Rapid approach to determine incompatibilities
- ▶ Multiple prototype formulations screened simultaneously
- ▶ Provides greater understanding of drug product stability
 - ▶ Impurity profile
- ▶ Allows rapid entry in full formulation development
- ▶ Predicts impact of container/closure on product stability

ASAP: Excipient Compatibility Study

- ▶ Sample count depends on study
 - ▶ 21 samples for 3 batches tested
- ▶ Amebis Incubation System
 - ▶ Incubation flask controls temperature/humidity



Case Study

- ▶ aNDA program
- ▶ Presence of lactose in reference listed drug (RLD)
 - ▶ Maillard reaction
 - ▶ Chemical reaction between amine and sugar
- ▶ One known degradation product
- ▶ Specification limit defined by ICH Q3B(R2) and monograph
- ▶ Drug product: oral tablet
- ▶ Strength: 0.5mg, 1mg, and 2mg
- ▶ Commercial packaging configuration: 100 and 1000 count HDPE bottle
- ▶ Timeline
 - ▶ 12 months to develop final formulation

Case Study

- ▶ API determined to be stable
- ▶ Three prototype formulations pressed into tablets
- ▶ Manufactured at bench scale
- ▶ Formulation 1 approximates RLD formulation

Prototype Formulation	1	2	3
Ingredient			
API	X	X	X
Microcrystalline Cellulose	X	X	X
Lactose	X		
Dicalcium Phosphate		X	
Starch 1500			X
Magnesium Stearate	X	X	X

Case Study

- ▶ ASAP design parameters

- ▶ 7 samples per prototype formulation
- ▶ Container closure system selected based on appropriate sized bottle for 100 tablets

Product Type:	Tablet
Length:	14 days
Number of Analyses:	7
Specification Limit:	0.2% label claim
Container:	75cc HDPE bottle/HIS cap

Case Study

- ▶ Testing design
 - ▶ Screening design for excipient selection
 - ▶ Minimal approach taken comparing three prototypes
 - ▶ Once prototype formulation selected and optimized then additional more robust ASAP study conducted


T (°C)	%RH	Time(Days)
50	75	14
60	50	14
70	10	14
70	75	1, 3
80	30	2

Case Study

- ▶ Stability indicating UHPLC method
- ▶ Reporting impurities equal or greater than 0.02% label claim
- ▶ 2 days for analytical analyses with a developed method

Mobile Phase	A:	0.1% TFA in H ₂ O	
	B:	Acetonitrile	
Gradient Program (Linear)	Time (min)	%A	%B
	Initial	90	10
	0.50	90	10
	5.00	25	75
	5.50	25	75
	5.55	90	10
	6.50	90	10
Column	Type:	Waters Acquity BEH C18, 50 x 2.1 mm, 1.7 μm particle size	
	Temperature:	40 °C	

Chromatogram



Raw Data: % Impurity

Formulation 1: Lactose

Deg Product	Control	50°C/75%RH 14 days	60°C/50%RH 14 days	70°C/11%RH 14 days	70°C/75%RH 1 days	70°C/75%RH 3 days	80°C/30%RH 2 days
Known 1	ND	0.032	0.078	0.076	0.039	0.065	0.099
RRT 3.21	ND	0.141	0.262	0.108	0.150	0.325	0.192
RRT 3.31	ND	0.058	0.256	0.064	0.086	0.059	0.179

Formulation 2: Dicalcium Phosphate

Deg Product	Control	50°C/75%RH 14 days	60°C/50%RH 14 days	70°C/11%RH 14 days	70°C/75%RH 1 days	70°C/75%RH 3 days	80°C/30%RH 2 days
Known 1	0.024	0.392	0.621	0.622	0.314	0.614	0.629

Formulation 3: Starch 1500

Deg Product	Control	50°C/75%RH 14 days	60°C/50%RH 14 days	70°C/11%RH 14 days	70°C/75%RH 1 days	70°C/75%RH 3 days	80°C/30%RH 2 days
Known 1	ND	0.055	0.036	0.043	0.046	0.115	0.040


Case Study

- ▶ Predicted results from ASAP®Prime meeting acceptance criteria of 0.2% label claim
 - ▶ All degradation products
- ▶ Packaging assessed
 - ▶ With and without desiccant
- ▶ ICH accelerated and CRT
- ▶ RRT 3.31 poor correlation
 - ▶ Potential secondary degradation


Formulation	1	1	1	2	3
Degradation Product	Known	RRT	RRT	Known	Known
Model correlation	0.96	0.97	0.62	0.99	0.95
Condition	% Probability of Passing				
40°C/75%RH 6 month no desiccant	98	90	77	0	99
40°C/75%RH 6 month 1g desiccant	99	100	90	36	100
25°C/60%RH 24 month no desiccant	99	100	90	73	99
25°C/60%RH 24 month 1 g desiccant	99	100	94	97	100
25°C/60%RH 36 month no desiccant	98	100	79	30	98
25°C/60%RH 36 month 1 g desiccant	99	100	87	80	99

Case Study

Formulation 2-Known-1 40° C/75% RH
no desiccant




Formulation 2-Known-1 25° C/60%RH
no desiccant



Case Study

- ▶ Formulation 1- RRT 3.31 40° C/75% RH
 - ▶ Higher error propagation



Jerusha Thompson

Case Study

- ▶ Results
 - ▶ Formulation 1 (RLD)
 - ▶ Two unique unidentified degradation products observed
 - ▶ Formulation 2
 - ▶ Known 1 would have difficulty meeting a shelf-life of 24 months without a desiccant
 - ▶ Failure at accelerated conditions would require additional testing at an intermediate condition
 - ▶ Formulation 3
 - ▶ Highest probability of passing all conditions without desiccant
- ▶ Study results will help justify the excipient and packaging selection in CMC section of submission

Alternative Analytical Technique

- ▶ Rapid entry into full formulation development with ASAP allows decreased timeline
 - ▶ Traditional approach to excipient studies ~3-6 months
 - ▶ ASAP ~ 2-4 weeks
- ▶ Rate limiting step can become analytical technology
 - ▶ Stability indicating HPLC method can take time to develop
 - ▶ ~1-2 months
 - ▶ Response of degradation products to parent compound critical
- ▶ Quantitative NMR (qNMR)
 - ▶ One method for all drug products
 - ▶ No method development required
 - ▶ No response of degradation product needed
 - ▶ No quantitative reference standard required


Alternative Analytical Technique

► qNMR Parameters

Field Strength	≥ 400-600 MHz
Nucleus	¹ H
Temperature	Regulated, typically 25-30°C
Number of Scans	≥ 32
Flip Angle	≤ 45°
Spectral Width (sw)	≥ 16 ppm
Relation Delay (D1)	30 s
Solvent	Chloroform
Acquisition Time	≥ 3 s

1. Pauli, G. F. et al. J. Med Chem. 2014, 57, 9220-9231
2. Weber, M. et al. J. Pharm. Biomed. Anal. 2014, 93, 102-110

Alternative Analytical Technique



Alternative Analytical Technique

- ▶ Predicted results for meeting acceptance criteria of 0.2%

Formulation	1	2	3
Degradation Product	Known 1	Known 1	Known 1
Model correlation	0.36	0.97	0.85
Condition	% Probability of Passing		
40°C/75%RH 6 month no desiccant	56.4	0.0	97.6
40°C/75%RH 6 month 1g desiccant	93.0	0.03	98.7
25°C/60%RH 24 month no desiccant	94.1	7.4	94.4
25°C/60%RH 24 month 1 g desiccant	98.8	60.2	96.0
25°C/60%RH 36 month no desiccant	83.5	0.48	91.0
25°C/60%RH 36 month 1 g desiccant	95.0	21.7	93.8

Alternative Analytical Technique

- ▶ Known 1 data for each prototype formulation analyzed using qNMR
- ▶ The results of the analysis predict prototype formulation 3 (Starch 1500) would have the greatest probability of meeting an acceptance criteria of 0.2% label claim
- ▶ Analysis of unknown degradation products not performed during this analysis
- ▶ Further optimization of extraction process necessary
- ▶ Both analytical techniques selected prototype formulation 3 as the most stable formulation

Conclusion

- ▶ The data generated from the excipient compatibility study predicted prototype formulation 3 would be the best candidate to move forward into full formulation development
- ▶ Analysis of each of the prototype formulation in the ASAP study also provided an early profile of likely degradation products present in each formulation
- ▶ Both analytical techniques investigating selected prototype formulation 3
 - ▶ qNMR analysis may allow for faster determination of instabilities in a drug product
- ▶ The predictions from the ASAP study will be further monitored through the project, and confirmed through traditional ICH stability studies
- ▶ ASAP approach to excipient compatibility testing allows faster entry into full formulation development ~2-5 month

Acknowledgement

- ▶ Patrick Nelson, Chemist II
- ▶ Limin Shi, Scientist II
- ▶ Sze Leung, Associate Chemist I
- ▶ Tyler Blanke, Associate Chemist I
- ▶ Byrant Gay, Ph.D., Scientist I
- ▶ Mark Gherke, Ph.D., Associate Director